

IMAPlate™ 5RC96

a Product from

Patented **I**ntelligent **M**ultifunctional **A**nalytical Technology

April 2012



-
- Introduction
 - Working principle
 - Classical liquid handling
 - Technical data
 - Application
 - Liquid transfer
 - Measurement
 - In-well homogeneous assay
 - In-well heterogeneous assay
 - 3D cell culture
 - Summary

A standard 96-well plate formatted
disposable miniature multi-utility lab device


with Integrated functions of

Pipette, Test Tube, Cuvette ...

The world's first analytical platform
capable of manually performing parallel, high throughput
liquid transfer, reaction, analysis and assay.





IMAPlate™ 5RC96 Working Principle – Capillarity

illustrations of unique self-dosed capillarity liquid handling (single unit)





Illustrations of Pipetting Solution to IMAPlate™ 5RC96

methods for pipetting 1- 5µl or <1µl of solution



methods for pipetting 15-25µl of solution



comparison between IMAPlate™ 5RC96 and Nunc. MaxiSorp F96 plate

	IMAPlate™ 5RC96	Nunc. MaxiSorp
Plate Format	SBS 96-well format	SBS 96-well format
Well Dimension (Φ x H)	1.1mm x 5mm	6.9mm x 10.7mm
Total volume	5µl (or 50µl)	360µl
Working volume (V)	5µl (or 15 – 25µl)	100µl
Surface area of working volume (Sa)	18mm ²	(60 + 34)mm ²
Surface area to volume ratio (Sa/V)	3.6mm ² /µl	0.94mm ² /µl
Initial binding velocity d[AB]/dt (B ⇔ Sa/V)	3.8 k _{on} [A][B]	1 k _{on} [A][B]
Complex density [AB] (B ⇔ Sa/V)	3.8 K _{eq} [A][B]	1 K _{eq} [A][B]
Diffusion distance from center to wall (X)	0.55mm	3.5mm
Diffusion time (X ² = D _{dif} •t)	0.3/D _{dif}	12/D _{dif}
Light path length	5 or 7.5mm (5 or 25µl)	3mm (100µl)

-
- ❑ 96 channel pipette & replicator for simultaneous liquid transfer
 - sampling, dilution, spotting, inoculation, ...

 - ❑ 96 micro-cuvette for UV-Vis-IR full spectrum measurement with easy sample recovery option
 - DNA, RNA and protein quantification at UV, compound analysis at UV-Vis-IR, assay result readout to enhance sensitivity.


 - ❑ 96 micro-well plate for reactions and assays
 - In-well homogeneous assays (e.g. miniature, high sensitivity assays)
 - In-well heterogeneous assays (e.g. miniature ELISA)
 - Solid phase enzymatic or binding assays for compound, ligand screen
 - 3D cell culture and cell based assays
 - Molecular biological assays
 - Proteomics applications
 - Multi-step solid phase chemical synthesis

-
- ❑ 96 channel pipette & replicator for simultaneous liquid transfer
 - sampling, dilution, spotting, inoculation, ...

 - ❑ 96 micro-cuvette for UV-Vis-IR full spectrum measurement with easy sample recovery option
 - DNA, RNA and protein quantification at UV, compound analysis at UV-Vis-IR, assay result readout to enhance sensitivity.

 - ❑ 96 micro-well plate for reactions and assays
 - In-well homogeneous assays (e.g. miniature, high sensitivity assays)
 - In-well heterogeneous assays (e.g. miniature ELISA)
 - Solid phase enzymatic or binding assays for compound, ligand screen
 - 3D cell culture and cell based assays
 - Molecular biological assays
 - Proteomics applications
 - Multi-step solid phase chemical synthesis

Illustrations of sampling, dilution, spotting, inoculation...






-
- ❑ 96 channel pipette & replicator for simultaneous liquid transfer
 - sampling, dilution, spotting, inoculation, ...

 - ❑ 96 micro-cuvette for UV-Vis-IR full spectrum measurement with easy sample recovery option
 - DNA, RNA and protein quantification at UV, compound analysis at UV-Vis-IR, assay result readout to enhance sensitivity.

 - ❑ 96 micro-well plate for reactions and assays
 - In-well homogeneous assays (e.g. miniature, high sensitivity assays)
 - In-well heterogeneous assays (e.g. miniature ELISA)
 - Solid phase enzymatic or binding assays for compound, ligand screen
 - 3D cell culture and cell based assays
 - Molecular biological assays
 - Proteomics applications
 - Multi-step solid phase chemical synthesis

Illustrations of assay in 96-well plate and measurement in IMAPlate™ 5RC96




perform assay in 96-well palte

transfer assay solution to IMAPlate

measure IMAPlate in a reader


The use of two wavelength absorbance measurement to calculate the absorbance for IMAPlate™ 5RC96



96-well plate measurement

$$\text{Absorbance} = -\log T = -\log(I/I_0) \quad (1)$$

The absorbance of the solution in 96-well plate is calculated according to the equation (1), where I_0 is the intensity of incident light and I is the intensity of transmitted light.



IMAPlate measurement

Due to small diameter of the capillarity chamber, the observed Abs is the sum of absorbance of the solution and blockage of incident light, which causes the baseline shifting upward. Assuming the blockage is the same over the spectrum, the absorbance can then be calculate by following equation (2).


$$\text{Absorbance} = \text{Abs}_{\text{Peak}} - \text{Abs}_{\text{Baseline}} \quad (2)$$

High sensitivity for absorbance and fluorescence measurement resulted from long light path-length of IMAPlate™ 5RC96

- ❖ Absorbance of **5µl** dye solution in IMAPlate™ 5RC96 = **60µl** in 384-well plate = **175µl** in 96-well plate
- ❖ Fluorescence count of **5µl** fluorescein solution in IMAPlate™ 5RC96 = **85µl** in 96-well plate


Absorbance vs Volume

(comparison of IMAPlate, 384- and 96-well plates)



Fluorescence vs Volume


(comparison of IMAPlate and 96-well plates)




* Data obtained from PerkinElmer EnSpire® Multimode Plate Reader.

High-throughput small sample volume quantification by touch-loading from a 96-well plate

DNA and protein measurement at UV





DNA (µg/ml)	Abs	SD	CV%
0	0.000	0.001	
3.6	0.028	0.003	10.9
7.3	0.071	0.003	4.3
14.5	0.160	0.001	0.6
29	0.329	0.004	1.2
58	0.632	0.006	1.0
116	1.239	0.015	1.2
232	2.449	0.004	0.2
23	0.256	0.011	4.3




BSA (mg/ml)	Abs	SD	CV%
0	0.000	0.001	
0.16	0.045	0.001	1.6
0.31	0.105	0.004	3.4
0.63	0.194	0.006	2.9
1.25	0.382	0.017	4.4
2.5	0.767	0.041	5.4
5	1.429	0.079	5.5
10	2.462	0.003	0.1
2	0.641	0.008	1.3

Good correlation with other measurement method

- DNA measurement by pipetting to IMAPlate
 - ❖ linear increase in Absorbance vs volume in the range of 1 to 5.5µl.
 - ❖ good correlation with nanodrop.



To increase sensitivity of ELISA by reducing substrate volume and transferring the solution to IMAPlate™ 5RC96 for measurement




$$\text{Abs} = \epsilon l c$$

- According to Beer-Lambert Law, light absorbance is the multiplication of absorption coefficient (ϵ), path length (l) and concentration (c).
- With same amount of enzyme in the well, the less volume of substrate solution, the higher concentration of colored product and the shorter path length for colored solution.
- Abs will not increase when using less substrate and measuring in the same well because the shorter path length cancels out the higher concentration of colored product.
- Once the colored solution is transferred to IMAPlate for measurement, the Abs increases due to the higher concentration of colored product and the longer path-length of IMAPlate.

Human IL6 ELISA


(15 μ l TMB & 5 μ l stop solution with touch-loading)




Spies P. et al. [A simple approach to improve the sensitivity of enzyme-linked immunosorbent assay: Using the IMAPlate 5RC96 for result readout](#), *Analytical Biochemistry*, 397 (2010): 48-55.

To increase sensitivity for commercial ELISA kit (pre-coated plate)

LEGEND MAX™ Human IL6 ELISA



LEGEND MAX™ Human IL17A ELISA



Procedure

- Prepare a set of extended standard with additional several lower concentration points.
 - Perform the ELISA according to the protocol in the kit **till to the step before the addition of substrate solution**. Then follow the procedure below.
- 1) *Add 40µl of TMB solution instead of 100µl to the wells of ELISA kit plate.*
 - 2) *Shake the plate on the plate shaker with 700-800 rpm for 15-20min.*
 - 3) *Pipette 10µl of the stop solution to the wells to stop the reaction.*
 - 4) *Transfer 40µl of the reaction mixture to the non-capillarity reservoir of the IMAPlate.*
 - 5) *Tap the IMAPlate several times horizontally to facilitate solution go to the capillarity chamber.*
 - 6) *Measure the IMAPlate at 450nm and 650nm in a plate reader.*
 - 7) *Calculate the true Abs by substrating Abs450 by Abs650 for plotting.*


-
- ❑ 96 channel pipette & replicator for simultaneous liquid transfer
 - sampling, dilution, spotting, inoculation, ...

 - ❑ 96 micro-cuvette for UV-Vis-IR full spectrum measurement with easy sample recovery option
 - DNA, RNA and protein quantification at UV, compound analysis at UV-Vis-IR, assay result readout to enhance sensitivity.


 - ❑ 96 micro-well plate for reactions and assays
 - In-well homogeneous assays (e.g. miniature, high sensitivity assays)
 - In-well heterogeneous assays (e.g. miniature ELISA)
 - Solid phase enzymatic or binding assays for compound, ligand screen
 - 3D cell culture and cell based assays
 - Molecular biological assays
 - Proteomics applications
 - Multi-step solid phase chemical synthesis

Illustrations of in-well homogeneous assay

5 μ l low volume assay with add-mix-measure protocol




15-25 μ l macro volume assay with add-mix-measure protocol with incubation



In-well miniature protein quantification assay


5µl Bradford assay

4µl of protein and 1µl of Bradford reagent mixed in the capillarity chamber by inverting the IMAPlate and measured immediately.



25µl BCA assay

2.5µl of protein and 22.5µl BCA working reagent incubated 30min at 37°C on a shaker in the non-capillarity reservoir of a parafilm sealed IMAPlate. Peeled off the parafilm and measured.




In-well miniature enzyme assay (15-25 μ l)



- e.g. *Candida Antarctica* Lipase B assay in IMAPlate for rapid colony screening

*M. A. Sciotti et. al. IMAPlate Based Miniature, High Sensitive, Rapid Screening Method for Detecting Bioengineered, Secreted Lipase Activities in Yeast Expression Systems. *CHIMIA*, Volume 64, No. 11, November 2010, Pages 789-792


- ✓ only 1/10 assay volume of 96-well plate
- ✓ higher sensitivity
- ✓ reducing assay cost
- ✓ minimizing upstream processing scale
- ✓ saving culture medium and reagents
- ✓ short screening cycle time




Illustrations of in-well miniature heterogeneous assay – ELISA



method for adding reagents and sample



method for washing




method for stopping reaction


large surface area to volume ratio and short diffusion distance leading to the increase in initial binding velocity and short time to reach equilibrium

- ❖ protein binding shown to obey single exponential kinetics
- ❖ formation of antigen antibody complex shown single exponential behavior as well


Time Course of



- 1) Coating: IgG-HRP in 1 $\mu\text{g/ml}$ BSA (diluted from 15 mg/ml BSA containing IgG-HRP)
- 2) TMB reaction: 5min for IMAPlate and 15 min for Nunc plate



- 1) Coating: 10 ng/ml GxR-IgG
- 2) Reaction: 1:20K R-IgG-HRP, 30min
- 3) TMB reaction: 15 min




- 1) Coating: 2 $\mu\text{g/ml}$ GxR-IgG, 30 min
- 2) Reaction: 1:150K R-IgG-HRP
- 3) TMB reaction: 5min for IMAPlate and 15 min for Nunc plate

In-well miniature ELISA


- ❖ small sample volume
- ❖ higher sensitivity
- ❖ short time to result
- ❖ less consumption of reagents

Rapid Miniature Troponin I ELISA
 IMAPlate: <3hrs incl. coating
 96-well plate: 7hrs + **ON** coating




from PerkinElmer Application Note: [Using the EnSpire Multimode Plate Reader to Measure IMAPlate-Based Rapid Miniature ELISA for the Quantification of Troponin I](#)

Human IL-6 ELISA



DuoSet kit from R&D systems: 1) coating O.N. 2) 1st incubation 120min 3) 2nd incubation 120min 4) SA-HRP 20min and 5) TMB 20min.

Human IL-17A ELISA



MAX™ Deluxe Kit from Biolegend: 1) coating O.N. 2) 1st incubation 120min 3) 2nd incubation 60min 4) SA-HRP 30min and 5) TMB 15min


Illustrations of in-well miniature heterogeneous assay – assays for compound screen

Compound screen with

solid phase enzyme assay

captive enzyme assay

competitive binding assay



Advantages


- ☐ easy wash procedure
- ☐ solid phase enzyme assay
 - able to eliminate the interference of colored compound on assay measurement
 - possible to reuse the same coated plate several times when enzyme is stable
- ☐ captive enzyme assay
 - reducing other enzymes interference
 - minimizing extract caused matrix effect for the assay

Why switching from 2D cell culture to 3D cell culture?



Drug candidates showing promise in monolayer 2D cell cultures often fail in later stages of drug development.

- cells altering the in vivo morphological and physiological characteristics when cultured in 2D
- constrained cell signaling mechanisms between adjacent cells due to limited cell-cell interactions in 2D
- 3D spheroids having similar cell-cell interactions like tissues in vivo
- 3D spheroid displaying very similar gene expression pattern of clinical specimens


Illustration of hanging drop 3 cell culture in Petri dish (limitations: medium change and spheroid transfer)




Illustrations of the preparation of spheroid micro-tissue and embryoid body in IMAPlate™ 5RC96



U87 spheroid in IMAPlate™ 5RC96
(2 days after seeding 5K, 2.5K and 1K cells: right to left)



LN319 spheroid in assay plate
(2 days after exposure to compound)



HEK Spheroid
(expressing GFP)

Advantages of using IMAPlate™ 5RC96 as a “hanging drop” cell culture device for 3D cell culture

1. easy to change cell culture medium
2. no risk of loss spheroid when aspirating the medium from the non-capillarity reservoir
3. rapid spheroid transfer from IMAPlate to 96- & 384-well plate
4. no harsh operating procedure for spheroid transfer
5. able to prepare and test spheroids in the same IMAPlate
6. capable of preparing large quantity of spheroids
7. easy to estimate the size of spheroids
8. cells readily to form spheroids in small capillarity chamber
9. possible to use a shaker to facilitate the formation of spheroids

Features of IMAPlate

- ❑ **Miniature:**
 - ✓ minimizing sample consumption
 - ✓ reducing assay cost
- ❑ **Unique and flexible liquid handling**
 - ✓ simultaneous filling/emptying samples
 - ✓ retaining pipette liquid handling options
 - ✓ avoiding air gap formation with open end
- ❑ **Improvement for reaction and analysis**
 - ✓ increasing sensitivity for measurement
 - ✓ rapid adsorption and binding
 - ✓ short time for reaching equilibrium
 - ✓ UV-Vis-IR full spectrum measurement
- ❑ **Compatible with 96-well plate readers and automated liquid handling systems**
 - ✓ no special equipments required

Benefits of using IMAPlate

- ✓ able to obtain multiparameter from limited amount of sample
- ✓ minimizing the scale of upstream process such as cell culture, fermentation, and protein purification
- ✓ no killing for small animal study and using less animal
- ✓ able to follow the same animal at time points and reducing variation of individual difference
- ✓ improving public images on animal research and industrial waste production (green technology)
- ✓ manual high-throughput
- ✓ short time to result
- ✓ reducing time for screening cycle
- ✓ enhancing lab productivity

NCL New Concept Lab GmbH

Eichenstrasse 22, CH4313 Moehlin, Switzerland

Tel: +41 61 853 08 20 / Fax: +41 61 853 08 23

Web: www.nclnewconceptlab.com



accompany you into the future